

# Rocky Mountain Conference on Magnetic Resonance

---

Volume 51 *51st Rocky Mountain Conference on Analytical Chemistry*

Article 1

---

7-19-2009

## 51st Rocky Mountain Conference on Analytical Chemistry

Follow this and additional works at: <https://digitalcommons.du.edu/rockychem>

 Part of the [Chemistry Commons](#), [Materials Science and Engineering Commons](#), and the [Physics Commons](#)

---

### Recommended Citation

(2009) "51st Rocky Mountain Conference on Analytical Chemistry," *Rocky Mountain Conference on Magnetic Resonance*: Vol. 51, Article 1.

DOI

<https://doi.org/10.56902/RMCMR.2009.51.1>

Available at: <https://digitalcommons.du.edu/rockychem/vol51/iss1/1>



This work is licensed under a [Creative Commons Attribution 4.0 International License](#).

This Conference Proceeding is brought to you for free and open access by Digital Commons @ DU. It has been accepted for inclusion in Rocky Mountain Conference on Magnetic Resonance by an authorized editor of Digital Commons @ DU. For more information, please contact [jennifer.cox@du.edu](mailto:jennifer.cox@du.edu), [dig-commons@du.edu](mailto:dig-commons@du.edu).

---

## 51st Rocky Mountain Conference on Analytical Chemistry

### Abstract

Final program, abstracts, and information about the 51st annual meeting of the Rocky Mountain Conference on Analytical Chemistry, co-endorsed by the Colorado Section of the American Chemical Society and the Society for Applied Spectroscopy. Held in Snowmass, Colorado, July 19-23, 2009.

### Copyright Statement / License for Reuse



This work is licensed under a [Creative Commons Attribution 4.0 International License](https://creativecommons.org/licenses/by/4.0/).

### Publication Statement

Copyright is held by the Rocky Mountain Conference on Magnetic Resonance. User is responsible for all copyright compliance.



## Unsurpassed Sensitivity...



## ...Enabling New Research Applications

Bruker BioSpin launches the world's first 263 GHz Solid State DNP-NMR spectrometer enabling DNP experiments at high field. The unique microwave source, an easy to use software-controlled high power gyrotron, has been specifically designed to provide you with a turn-key solution for DNP-enhanced solids NMR experiments. Optimum beam propagation to the sample is ensured by microwave transmission lines which results in maximum polarization efficiency. Together with our new low-temperature MAS probe technology, the 263 GHz DNP accessory delivers unsurpassed sensitivity enabling exciting new applications for solid-state NMR researchers.

[www.bruker-biospin.com/dnp-solids](http://www.bruker-biospin.com/dnp-solids)

think forward


Bruker BioSpin pioneered the world's first commercial mm-wave 263 GHz EPR spectrometer, ELESYS E780, representing a first step for Bruker's EPR division into quasi-optical microwave technology. It incorporates a unique superconducting magnet that can be ramped up to 12T and is combined with new probe technology for optimum sensitivity, even on large samples up to 5 mm. Based on the well-proven Bruker ELESYS concept it provides multiple turn-key operation modes including, CW-, Pulse-EPR, ENDOR- and ELDOR, thus enabling research groups for the first time to routinely use Very High Frequency EPR technology.

[www.bruker-biospin.com/e780](http://www.bruker-biospin.com/e780)

51<sup>ST</sup> ROCKY MOUNTAIN CONFERENCE ON ANALYTICAL CHEMISTRY

July 19 – 23, 2009

# 51<sup>ST</sup> ROCKY MOUNTAIN CONFERENCE ON ANALYTICAL CHEMISTRY



## FINAL PROGRAM AND ABSTRACTS

ENDORSED BY:

Colorado Section – American Chemical Society  
& Society for Applied Spectroscopy

July 19 – 23, 2009

Snowmass Conference Center • Snowmass, Colorado, USA

[www.rockychem.com](http://www.rockychem.com)